

Can DNA-metabarcoding of multi-species samples from routine egg surveys be used to monitor the ichthyoplankton assemblage?


Associate Professor Tim Ward (IMAS)

S. Appleyard, S. Catalano, F. Bailleul, A. Ivey,
A. McGrath, A. Miskiewicz, G. Grammer

SPF May 2026




Broad-scale Egg Surveys



- ~45 surveys of SE Australia since 1995
- Each survey covers 60,000 to 150,000 km² and includes hundreds of samples
- sardine, jack mackerel, redbait, blue mackerel
- Used to set TACs in SA Sardine Fishery & Commonwealth SPF (TACs~80,000 t)
- **Can we value-add?**

Can we use DNA metabarcoding to identify fish species in a plankton sample?

- Samples from 4 surveys
 - different handling methods
 - different species assemblages (seasons & depths)
- Fish larvae and eggs separated
- Morphological identification
- Tested 2 mtDNA-metabarcoding assays
 - Assays target different ends of COI gene
 - Cytochrome oxidase sub-unit I (COI)
 - COI_E (Shokralla et al. 2015):
 - LH (5') end
 - COI_L (Leray et al. 2013):
 - RH (3') end




Bulked Larval Libraries

Detailed morpho


R41-2-L

Species	count	Mean size (mm)	characteristics
A	11	9	Gonostomatidae
B	4	12.5	Myctophidae sp. 1
C	5	10	Myctophidae sp. 2
D	7	8.5	Myctophidae sp. 3
E	2	6	Myctophidae sp. 4
F	1	8	Myctophidae sp. 5
G	1	11	Myctophidae sp. 6
H	1	7	Myctophidae sp. 7
I	1	7	Carapidae? – can't see rostrum behind head tho
J	1	4	? Playcephalid?
K	2	4.5	Sebastidae – Helicolenus sp.
L	1	4	?
M	1	3	Yolk sac?
N	1	4.5	Stargazer?
O	2	6.5	Myctophidae sp. 8
P	1	6	Myctophidae sp. 9
Q	2	8	Myctophidae sp. 10

N=44



DNA extraction



'Rules' developed iteratively for filtering and acceptance of presence

- 98-100% pairwise identity
- >20 reads unless

R44-1-L

Species	count	Mean size (mm)	characteristics
A	1	9	Myctophidae


N=1



Shokralla



Leray




- Different numbers of individuals bulked in sample

Bulked Egg Libraries


Detailed morpho

R5-2-E




Species	count	Mean size (µm)	characteristics
A	607	0.9	Many species - very diverse. Possibly Tracharia, early etc, plus several others. All between 0.8-1.0µm

N=637



R23-2-E



Species	count	Mean size (µm)	characteristics
A	1	1.4	
B	1	1.0	
C	1	1.0	

N=3



Shokralla




Leray

DNA extraction


- Different numbers of individuals bulked in samples

Comparisons: Morphology vs Leray & Shokralla

Eggs



Larvae



	Family	Genus	Species
Morpho	6	6	6
Leray	84	154	203
Shokralla	89	167	222
TOTAL (L + S)	91	177	248

Assays detected ~200 species > morphology

	Family	Genus	Species
Morpho	77	52	29
Leray	87	178	236
Shokralla	88	188	250
TOTAL (L + S)	90	198	275

Why?

Identified by MORPHO but not by either Assay:

- 8 Families
- 10 Genera
- 8 Species



All Families, Genera & Species identified by MORPHO were ID'd by both assays

Checking Comparisons: Larvae Morphology vs Leray & Shokralla

8 Families identified by MORPHO but not by Assays

Family	Morpho	Leray	Shokralla
Carapidae	Y	-	-
Champsodontidae	Y	-	-
Chlorophthalmidae	Y	-	-
Howellidae	Y	-	-
Ipnopidae	Y	-	-
Lophiidae	Y	-	-
Nemichthyidae	Y ?	-	-
Syngnathidae	Y	-	-

10 Genera identified by MORPHO but not by Assays

Genus	Morpho	Leray	Shokralla
Bathylagus	Y	-	-
Champsodon	Y	-	-
Chlorophthalmus	Y	-	-
Gymnapistes	Y	-	-
Hime	Y	-	-
Howella	Y	-	-
Melanostomias	Y	-	-
Neosebastes	Y	-	-
Rhombosolea	Y	-	-
Stigmatopora	Y	-	-




8 Species identified by MORPHO but not by Assays

Species	Morpho	Leray	Shokralla
<i>Aldrichetta fosteri</i>	Y ?	-	-
<i>Centroberyx affinis</i>	Y ?	-	-
<i>Cepola australis</i>	Y	-	-
<i>Cynoglossus broadhursti</i>	Y ?	-	-
<i>Gaidropsarus novaezealandiae</i>	Y ?	-	-
<i>Gymnapistes marmoratus</i>	Y	-	-
<i>Rhombosolea tapirina</i>	Y ?	-	-
<i>Trachurus novaezealandiae</i>	Y	-	-

Reasons:
Morpho: mis-ID's

Assays: <98% identity or low detections (<20)

Species record not in reference database




Checking Comparisons: Eggs -> Morphology vs Leray & Shokralla

All Families, Genera and Species ID'd by MORPHO ID'd by both Assays


Family	Genus	Species	Morpho	Leray	Shokralla
Bothidae			Y	Y	Y
Carangidae	<i>Trachurus</i>	<i>Trachurus declivis</i>	Y	Y	Y
Clupeidae	<i>Sardinops</i>	<i>Sardinops sagax</i>			
	<i>Etrumeus</i>	<i>Etrumeus teres</i>	Y	Y	Y
Emmelichthyidae	<i>Emmelichthys</i>	<i>Emmelichthys nitidus</i>	Y	Y	Y
Engraulidae	<i>Engraulis</i>	<i>Engraulis australis</i>	Y	Y	Y
Scombridae	<i>Scomber</i>	<i>Scomber australasicus</i>	Y	Y	Y

DNA good for eggs: majority of fish eggs really hard to identify morphologically and many undescribed


SA Sardine Survey



SA Sardine Survey




Blue Mackerel East Survey




Eggs

Comparisons: Leray vs Shokralla - % identification


Family



Genus




Species




Larvae


Family



Genus





Species



Comparisons: Morphology Counts vs Reads (Leray)

Sardine (*Sardinops sagax*) - Eggs



Reasons:


- Inhibition
- Swamping

Mismatch between morphological counts and number of reads demonstrates approach is not quantitative


Practical applications: Species assemblages)

Both assays useful for identifying biodiversity hotspots and potential monitoring changes over time.

Eggs




Larvae



Practical applications: Estimating Spawning Area for the DEPM and other specie

Sardinops sagax Sardine



A in 2020 Sardine DEPM

~82,000 km² (Ward et al. 2020)

Spawning area:

- Voronoi natural neighbour method
- Presence/absence of eggs in grid cell


A in 2005 Blue Mackerel DEPM

~36,000 km² (Ward & Rogers 2007)

Eggs

Scomber australasicus

Blue Mackerel



Key Findings

Morphological identification:

- Advantage: counts not presence/absence
- Disadvantages:
 - difficult, time consuming, expensive
 - can be prone to error
 - <50% identification of spp for larvae, much less spp ID for eggs

DNA-metabarcoding (Leray and Shokralla assays):

Advantages:

- 2 assays very similar results
 - Shokralla generally > number species detected,
 - BUT one assay may detect a target species better than other (e.g. Trachurus)
- more spp and genera identified compared to morphology,
- relatively cheaper
- does not require (rare) high level taxonomic expertise
- good for eggs – estimating spawning area
- Handling protocol (ultraclean vs standard) -- less important than initially thought
- **Hybrid approach (coarse morpho + DNA-mb) → Valuable**
- **Funding: despite surveys being cost-recovered, govt funding limited**

